

# Human Bladder Microvascular Endothelial Cell cDNA (HBdMEC cDNA (Formerly known as HBlMEC cDNA)) Catalog #4304

Volume: 20 µl (approximately 20 reactions)

# **Product Description** Human Bladder Microvascular Endothelial Cell PCR-ready first strand cDNA (HBdMEC cDNA) is prepared from RNA extracted from Human Bladder Microvascular Endothelial Cells by a Qiagen RNeasy kit. 2 $\mu$ g total RNA is then reverse-transcribed using an Applied Biosystems' High-Capacity cDNA Reverse Transcription kit. 1 $\mu$ l cDNA is sufficient for one PCR reaction. cDNA from ScienCell Research Laboratories is convenient and cost effective for researchers as it eliminates the need to acquire expensive tissues.

# **Quality Control**

- The  $A_{260}/A_{280}$  ratio of the cDNA is measured to be between 1.7-1.9 using a Beckman Coulter DU 800 spectrophotometer.
- The cDNA is performance-tested in a control PCR using  $\beta$ -actin primers.

## Application

- PCR
- DNA microarray and gene expression profiling
- Cloning of specific genes
- Gene alternative splicing and mutation detection

## **Product Use**

HBdMEC cDNA is for research use only. It is not approved for human or animal use, or for application in *in vitro* diagnostic procedures.

## Storage

Store cDNA at -20°C upon arrival; aliquot the cDNA to avoid repeated freeze/thaw cycles.

## Shipping

Dry ice