

# **Meningeal Cell Growth Supplement** (MCGS)

Catalog #1452

# **Product Description**

Meningeal Cell Growth Supplement (MCGS) is a sterile, concentrated (100X) solution which contains growth factors, hormones, and proteins necessary for the culture of normal meningeal cells. The supplement is formulated (quantitatively and qualitatively) to provide an optimally balanced environment that selectively supports the growth of normal meningeal cells *in vitro*.

#### **Product Use**

MCGS is for research use only. It is not approved for human or animal use, or for application in *in vitro* diagnostic procedures.

### Storage

Store MCGS at -20°C; once added to medium, store at 4°C and protect from light. Do not refreeze after thawing.

## Shipping

Dry ice.

#### **Instructions for Use**

For best results, use MCGS in conjunction with 500 ml Meningeal Cell Medium-b (MCM-b, Cat #1401-b) and 10 ml fetal bovine serum (FBS, Cat #0010).

Thaw MCGS at 37°C. Gently tilt the tube several times to ensure complete mixing. Spray the medium bottle and supplement tube with 70% ethanol and wipe to remove excess liquid. In a sterile field, remove the caps without touching the interior threads with fingers. Add MCGS to the basal medium and mix well. Since several components are light-labile, the medium should not be exposed to light for extended periods. We do not recommend warming medium in a 37°C water bath prior to use. When stored in the dark at 4°C, the reconstituted medium is stable for one month.

Caution: If handled improperly, some components of the medium may present a health hazard. Take appropriate precautions when handling it, including the wearing of protective clothing and eyewear. Dispose of properly.