



## Human Epidermal Keratinocyte-adult MicroRNA (HEK-a miRNA)

Catalog #2117, 2117-2, or 2117-5

Quantity: 1  $\mu$ g, 2  $\mu$ g, or 5  $\mu$ g

### Product Description

Human Epidermal Keratinocyte-adult microRNA (HEK-a miRNA) is prepared from early passage Human Epidermal Keratinocytes-adult using Life Technologies' *mirVana*<sup>™</sup> miRNA Isolation Kit. The microRNA is purified by organic extraction and enriched by immobilization of RNA on glass-fiber filters. The microRNA is eluted and stored in nuclease-free water. MicroRNA from ScienCell Research Laboratories is a convenient and cost effective alternative to acquiring expensive tissues.

### Quality Control

- The  $A_{260}/A_{280}$  ratio of the microRNA is measured to be 1.8-2.1 using Beckman Coulter DU 800 spectrophotometer.

### Application

- Solution hybridization assays
- Blot hybridization
- Microarray analysis
- Control miRNA profile

### Product Use

HEK-a miRNA is for research use only. It is not approved for human or animal use, or for application in *in vitro* diagnostic procedures.

### Storage

Store the microRNA immediately at  $-80^{\circ}\text{C}$  upon arrival; aliquot the microRNA to avoid repeated freeze/thaw cycles.

### Shipping

Dry ice.