



## **Human Non-Pigmented Ciliary Epithelial Cell MicroRNA (HNPCEpiC miRNA)**

Catalog #6587, 6587-2, or 6587-5

Quantity: 1 µg, 2 µg, or 5 µg

### **Product Description**

Human Non-Pigmented Ciliary Epithelial Cell microRNA (HNPCEpiC miRNA) is prepared from early passage Human Non-Pigmented Ciliary Epithelial Cells using Life Technologies' *mirVana*<sup>™</sup> miRNA Isolation Kit. The microRNA is purified by organic extraction and enriched by immobilization of RNA on glass-fiber filters. The microRNA is eluted and stored in nuclease-free water. MicroRNA from ScienCell Research Laboratories is a convenient and cost effective alternative to acquiring expensive tissues.

### **Quality Control**

- The  $A_{260}/A_{280}$  ratio of the microRNA is measured to be 1.8-2.1 using Beckman Coulter DU 800 spectrophotometer.

### **Application**

- Solution hybridization assays
- Blot hybridization
- Microarray analysis
- Control miRNA profile

### **Product Use**

HNPCEpiC miRNA is for research use only. It is not approved for human or animal use, or for application in *in vitro* diagnostic procedures.

### **Storage**

Store the microRNA immediately at -80°C upon arrival; aliquot the microRNA to avoid repeated freeze/thaw cycles.

### **Shipping**

Dry ice.