

Recombinant Murine basic Fibroblast Growth Factor (rmbFGF)

Catalog Number: 124-02

Description Basic Fibroblast Growth Factor (bFGF) is a member of the FGF family of mitogenic peptides

which is comprised of at least 22 proteins showing 35-55% amino acid sequence conservation. Unlike other FGF family members, bFGF and acidic FGF (aFGF) lack signal peptides and are secreted via a different mechanism other than the classical protein secretion pathway. The predicted 17 kDa bFGF isoform can be located in both the cytoplasm and the nucleus. The use of alternative transcription start sites within the gene produces several different isoforms

ranging from 21-24kDa proteins.

Synonyms FGF-2, Fibroblast Growth Factor-basic

AA Sequence MPALPEDGGAAFPPGHFKDPKRLYCKNGGFFLRIHPDGRVDGVREKSDPHVKLQLQA

EERGVVSIKGVCANRYLAMKEDGRLLASKCVTEECFFFERLESNNYNTYRSRKYSSW

YVALKRTGQYKLGSKTGPGQKAILFLPMSAKS

Source Escherichia coli

Molecular Weight Approximately 16.2 kDa, a single non-glycosylated polypeptide chain containing 146 amino

acids.

Purity >98% by SDS-PAGE and HPLC analyses.

Biological Activity Fully biologically active. The ED₅₀ is < 0.1ng/ml, corresponding to a specific activity of 1 x

10⁷ units/mg, as determined by proliferation of BaF3 cells expressing FGF receptors

Physical Appearance White lyophilized powder.

Formulation Lyophilized from a 0.2µm filtered solution in PBS, pH 7.4.

Endotoxin $< 1EU/\mu g$ of growth factor as determined by LAL method.

Reconstitution Reconstitute in sterile distilled water containing 0.1% BSA to a concentration of 0.1-1.0

mg/mL.

Storage Storage Store at -20°C after receiving. Upon reconstitution, store at 2-8°C for up to one week. For

maximal stability, aliquot and store at -20°C. Avoid repeated freeze/ thaw cycles.

Usage This product is for research use only. It is not approved for use in humans, animals, or *in vitro*

diagnostic procedures.